



Case Study #5: Imprecision of Repetitions at Levels 2 through 5

Initial Results: A laboratory performed routine calibration verification / linearity testing using VALIDATE[®] GC3. Results for all enzymes demonstrated imprecision at levels 2 through 5. The following AST report was generated using MSDRx[®], LGC Maine Standards Data Reduction software:

GC 3 AST

suggested total allowable error is 5 U/L or 20%, whichever is greater								Linearity Across the Reportable Range: GC 3 AST						
LX	Rep 1	Rep 2	Rep 3	Accept Comments						<u> </u>				
B N/A	A 1	2	2					1,200	ł:				X	
1 1.0) 4	5	4	Tested 4.33 to 1,057.7 U/L					1,000	4				
2 2.0	363	305	270	Validated to U/L				"	800	1				
3 3.0	0 675	602	529	Mean versus Target Regression y = 0.874x + 29.979					Value	1			-	
4 4.0	0 868	805	(80)							1:				
5 5.0	0 1095	1027		r^2=0.9843 SEy,x=50.1					400	+				
X	Target	Mean	+/- Diff	f % Diff	+/- Limit	% Limit			200	1				
1.0	7.50	4.33	** 3.1	7 42.3%	2.50	N/A			0		·			
2.0	306.33	312.67	6.3	4 2.1%	30.63	10%				1	2	3	4	5
3.0	605.17	602.00	3.1	7 0.5%	60.52	10%					Linea	ar Dilution	Level	
4.0	904.00	817.67	86.3	3 9.5%	90.40	10%			ſ	* Target	– Lim	it 🔳	Recovered	
5.0	1,202.83	1,057.67	145.1	6 ** 12.1%	120.28	10%			L					

The laboratory contacted LGC Maine Standards Technical Support. Technical Support advised the laboratory that the results were not consistent with Peers and that the curve does not demonstrate typical reproducibility of repetitions of levels 2 through 5.

Troubleshooting: The laboratory took the troubleshooting step of contacting the instrument manufacturer and requested a service call. During the service call, a cracked sample syringe case was discovered and replaced. To confirm the issue was corrected the laboratory re-ran their VALIDATE[®] GC3 test kit. The updated MSDRx[®] report for AST showed that all levels were within the statistical limits and precision of the repetitions at all levels acceptable. The laboratory accepted the updated results and determined that they had validated the linearity across the reportable range of the method.

GC 3 AST

suggested total allowable error is 5 U/L or 20%, whichever is greater							Linearity Acr	Linearity Across the Reportable Range: GC 3 AST					
L	Х	Rep 1	Rep 2	Rep 3		-+ □	C						
В	N/A	1	2	2		pt	Comme	1,000					
1	1.0	3	4	4	Tested 3.67 to 1,019.7 U/L			800					
2	2.0	247	248	247	Validated <u>3.67</u> to <u>1.019.7</u> U/L <u>Mean versus Target Regression</u> y = 1.035x - 5.043				1				
3	3.0	499	500	501							····		
4	4.0	784	786	782				S 400					
5	5.0	1016	1021	1022	r^2=0.9992 SEy,x=10.8								
X	Tar	rget	Mean	+/- Diff	% Diff	+/- Limit	% Limit	200	A Statement of the second s				
1.0) :	2.17	3.67	1.50	69.1%	2.50	N/A	0-1			· · · · · ·		
2.0) 25	0.33	247.33	3.00	1.2%	25.03	10%	1	2	3 4	5		
3.0) 49	8.50	500.00	1.50	0.3%	49.85	10%		Linear D	ilution Level			
4.0) 74	6.67	784.00	37.33	5.0%	74.67	10%	¥ Ta	arget – Limit	Recovered			
5.0	99	4.83 1	,019.67	24.84	2.5%	99.48	10%						

Summary: LGC Maine Standards Technical Support advised that recovery and reproducibility was not consistent with Peers. In this case, a cracked syringe case caused the impression seen between repetitions. This case study emphasizes that calibration and running QC may not always detect instrument issues.