

## **Malidate**

## Case Study #4: Low Bilirubin Recovery at All Levels and Results Appear Visually Nonlinear

**Initial Results:** A laboratory performed routine calibration verification / linearity testing using VALIDATE<sup>®</sup> GC 4. Results for both Total Bilirubin and Direct Bilirubin were well below the upper range of the method, and were visually nonlinear. The following TBIL report was generated using MSDRx<sup>®</sup>, the LGC Maine Standards Data Reduction software:

## GC 4 TBIL

published CLIA total allowable error is 0.4 mg/dL or 20%, whichever is greater

mg/dL or 20%, whichever is greater								Linearity A	Across th	e Reportable	Range: G	C 4 TBIL			
L	Х	Rep 1	Rep 2	Rep 3				monto	7						
В	N/A				יא בון	ccept			6						
1	1.0	0.2	0.2	0.0	Tested	0.133 t	o 6.73 mg	J/dL	5						
2	2.0	1.2	1.1	1.2	Validate	ed	to	mg/dL					1		
3	3.0	2.5	2.5	2.5		ersus Targ	ot Pograd		Values					5	
4	4.0	4.4	4.3	4.3		83x - 0.415		ssion	<b>R</b> 3						
5	5.0	6.8	6.7	6.7	r^2=0.9	721 SEy,x	=0.4		2						
X	Targ	get Me	ean +/-	Diff 9	6 Diff	+/- Limit	% Limit	]	1						
1.0	0.0	083 0.1	133 0	.050	60.2%	0.200	N/A	1							
2.0	1.2	267 1.1	167 0	.100	7.9%	0.200	N/A	1	0-	· · · ·	+ •	· · · ·		· · ·	
3.0	2.4	150 2.	500 0	.050	2.0%	0.245	10%	]	1		2 L	3 inear Dilutio	n Level	4	5
4.0	3.6	633 4.3	333 0	.700 **	19.3%	0.363	10%	]							
5.0	4.8	317 6.	733 1	.916 **	39.8%	0.482	10%	] [		L	<b>X</b> Target	- Limit	Recov	erea	

The laboratory contacted LGC Maine Standards Technical Support. Technical Support advised the laboratory that the results were not consistent with Peers and that the curve was typical of product that was not stored properly. The proper storage conditions for VALIDATE<sup>®</sup> GC 4 is -10 to -25°C in a non-frost-free freezer. More details regarding the need for non-frost-free freezers are provided under General Technical Service Bulletins located at www.mainestandards.com/support/tech-bulletins.php.

**Troubleshooting:** The laboratory took the troubleshooting step of reordering a fresh kit of GC 4 and rerunning the calibration verification / linearity experiment. The updated MSDRx<sup>®</sup> report for TBIL shows that all Levels are within the statistical limits. The laboratory accepted the updated results and determined that they had validated the linearity across the reportable range of the method.

GC 4 TBIL

published CLIA total allowable error is 0.4 mg/dL or 20%, whichever is greater

mg/d	L or 20	0%, which	never is g	reater				Linearity Across the Reportable Range: GC 4 TBIL
L	х	Rep 1	Rep 2	Rep 3		ccept		
В	N/A				יי טך		_	25
1	1.0	0.3	0.2	0.2	Tested	0.233 t	o 27.30 m	
2	2.0	6.4	6.4	6.4	Validate	ed 0.233 t	o 27.30 m	20
3	3.0	12.9	12.9	13.0			at Deere	15-
4	4.0	19.5	19.5	19.6		<u>ersus Targ</u> 59x - 0.356		- Cal
5	5.0	27.1	27.6	27.2		980 SEy,>		10
X	Targ	get Me	an +/-	Diff	% Diff	+/- Limit	% Limit	5
1.0	0.1	72 0.2	233 0	0.061	35.5%	0.200	N/A	
2.0	6.5	522 6.4	100 C	).122	1.9%	0.652	10%	
3.0	12.8	372 12.9	933 0	0.061	0.5%	1.287	10%	1 2 3 4 5 Linear Dilution Level
4.0	19.2	222 19.5	533 (	0.311	1.6%	1.922	10%	
5.0	25.5	572 27.3	300 1	.728	6.8%	2.557	10%	★ Target - Limit Recovered

Summary: LGC Maine Standards Technical Support advised that recovery was not consistent with Peers. In this case, improper storage damaged the VALIDATE<sup>®</sup> test set. This case study emphasizes the need to store products per the manufacturer's instructions. Storage requirements for VALIDATE<sup>®</sup> products are printed on the kit boxes, bottle labels and in the package inserts.