

TECHNICAL SERVICE BULLETIN



Technical Service Bulletin #04-003

Date: September 28, 2004
Re: VALIDATE[®] TDM 1
P/N: 301ri
Subject: Gentamicin / Tobramycin Recovery on Hitachi Systems

Ongoing quality assurance and customer feedback have shown that Gentamicin and Tobramycin recovery in VALIDATE[®] TDM 1 calibration verification / linearity test set, Product #301ri, are not compatible with Roche Catalog numbers 11815342 (Gentamicin II) and 11815385 (Tobramycin II) for the Hitachi instrument systems.

Based on this information, the use of this product on instrument systems utilizing these reagents for gentamicin and/or tobramycin is not recommended.

Please refer any questions related to this bulletin to our Technical Service department.

Thank you,

Technical Support Department
Maine Standards Company
765 Roosevelt Trail
Windham, ME 04062
1-207-892-1300
1-800-377-9684
Fax: 1-207-892-2266
E-mail: info@mainestandards.com

TSB #04-003