

TECHNICAL SERVICE BULLETIN



Technical Service Bulletin #05-003

Date: July 22, 2005
Re: VALIDATE[®] UC2
P/N: 702bc
Subject: Urine Calcium (CALC) Reagent Incompatibility

Ongoing quality assurance and customer feedback have shown that the VALIDATE[®] UC2 calibration verification / linearity test set, Reorder number 702bc, is incompatible with Beckman-Coulter **CALC reagent Part Numbers 450214 and 472095**. Assaying this product for Urine Calcium with these reagents will result in Suppressed Values for most Levels.

Urine Calcium in VALIDATE[®] UC2 is compatible with Beckman-Coulter reagent **Part Numbers CA 42755 and CA3 450890**.

Please refer any questions related to this bulletin to our Technical Support department.

Thank you,

Technical Support Department
Maine Standards Company
765 Roosevelt Trail
Windham, ME 04062
1-207-892-1300
1-800-377-9684
Fax: 1-207-892-2266
E-mail: info@mainestandards.com