

## TECHNICAL SERVICE BULLETIN



### Technical Service Bulletin #05-006

Date: December 8, 2005  
Re: VALIDATE<sup>®</sup> TDM 1  
P/N: 301ab  
Subject: Gentamicin / Tobramycin recovery

Ongoing quality assurance and customer feedback have shown that Gentamicin and Tobramycin recovery in VALIDATE<sup>®</sup> TDM 1 calibration verification / linearity test set, Order No. 301ab, are not compatible with the Microgenics methods.

Based on this information, the use of this product on instrument systems utilizing Microgenics reagents for gentamicin and/or tobramycin is not recommended.

Please refer any questions related to this bulletin to our Technical Support department.

Thank you,

Technical Support Department  
Maine Standards Company  
765 Roosevelt Trail  
Windham, ME 04062  
1-207-892-1300  
1-800-377-9684  
Fax: 1-207-892-2266  
E-mail: [info@mainestandards.com](mailto:info@mainestandards.com)