



Technical Service Bulletin #07-001

Date: January 15, 2007
Re: VALIDATE[®] Urine Chem 4 (UC4) bc
P/N: 704bc
Subject: Urine Calcium (CALC) Announcement
Reference: TSB #05-003, July 22, 2005

The Urine Chemistry 4 (UC4) bc product is now compatible with Beckman Coulter[®] CALC Reagent Part Numbers 450214 and 472905. Urine Calcium in VALIDATE[®] UC4bc also continues to be compatible with Beckman Coulter[®] reagent Part Numbers CA 442755 and CA3 450890.

The improvements will be implemented with Kit Lot # 74AA00307 Exp 2008-03-04 and newer.

Please refer any questions related to this bulletin to our Technical Support department.

Thank you,

Technical Support Department
Maine Standards Company
765 Roosevelt Trail
Windham, ME 04062
1-207-892-1300
1-800-377-9684
Fax: 1-207-892-2266
E-mail: techsupport@mainestandards.com