



Technical Service Bulletin #08-002

Date: March 3, 2008
Re: VALIDATE[®] Urine Chem 5 (UC5) bc
P/N: 705bc
Subject: Urine Amylase (AMY) Announcement

The Amylase in the VALIDATE[®] UC5 bc product is now designed to test the reportable range of the Beckman AMY reagent P/N 442775. The typical range of Amylase in the UC5 bc product is 5 – 800 U/L

The improvements will be implemented with Kit Lot # 75AM04908 Exp 2009-02-17 and newer.

Please refer any questions related to this bulletin to our Technical Support department.

Thank you,

Technical Support Department
Maine Standards Company
765 Roosevelt Trail
Windham, ME 04062
207-892-1300
800-377-9684
Fax: 207-892-2266
E-mail: techsupport@mainestandards.com