

**Date Issued:** September 30, 2009

**Product:**


<b>VALIDATE<sup>®</sup> Product</b>	<b>REF</b>
<b>TDM1vt</b>	<b>301vt</b>

**Subject:** Addition of Analytes to TDM1vt

The TDM1vt Product has been improved to test the reportable range of more TDMs on the Ortho-Clinical Vitros and Vitros 5,1 systems.

<b>TDM 1: Order No. 301vt, Five levels, 3 mL per level Storage: -10 to -25°C</b>	
<b>Analyte</b>	<b>Range*</b>
ACTM (Acetaminophen)	10 – 200 µg/mL
CARB (Carbamazepine)	3.0 – 20.0 µg/mL
DIGN (Digoxin)	0.40– 4.00 ng/mL
GENT (Gentamicin)	0.6 – 10.0 µg/mL
PHNO (Phenobarbital)	3.0 – 80.0 µg/mL
PHYT (Phenytoin)	3.00– 40.00µg/mL
SALY (Salicylate)	1.0 – 40.0 mg/dL
THEO (Theophylline)	1.00– 40.00µg/mL
TOB (Tobramycin)	0.6 – 10.0 µg/mL
VALP (Valproic Acid)	10 – 150.0µg/mL
VANC (Vancomycin)	5.0 – 50.0 µg/mL

The improvements will apply beginning with the lot number below.

<b>VALIDATE<sup>®</sup> Product</b>	<b>REF</b>	<b>LOT</b>	 <b>Expiration Date</b>
TDM1vt	301vt	30AB27109	2010-09-28

Please refer any questions related to this bulletin to our Technical Support department.

**Contact Information:**

Customer Support Department  
 Maine Standards Company  
 765 Roosevelt Trail  
 Windham, ME 04062  
 Phone: 1-800-377-9684  
 Fax: 1-207-892-2266  
 E-mail: [Info@mainestandards.com](mailto:Info@mainestandards.com)