




TECHNICAL SERVICE BULLETIN
#09-004

Date Issued: October 13, 2009

Product:

VALIDATE[®] Product	REF	LOT	 Extended Storage Temperature
Chem 20	200	All Lots	-10° to -25°C
Chem 21	201	All Lots	-10° to -25°C
CM2	402bc	All Lots	-10° to -25°C
TDM1	301ab	All Lots	-10° to -25°C
TDM1	301db	All Lots	-10° to -25°C
TDM1	301vt	All Lots	-10° to -25°C
LP	501bc	All Lots	-10° to -25°C
LP	501db	All Lots	-10° to -25°C
LP	501ro	All Lots	-10° to -25°C
GC1	1100db	All Lots	-10° to -25°C
GC4	1400db	All Lots	-10° to -25°C

Subject: Extended Storage Temperature for VALIDATE[®] test kits

Explanation: Ongoing Quality Assurance testing has led us to extend the storage temperature for all lot numbers of VALIDATE[®] calibration verification / linearity test kits listed above from -10° to -20°C to -10° to -25°C.

Please retain this bulletin for your records.

Contact Information:

Customer Support Department
Maine Standards Company
765 Roosevelt Trail
Windham, ME 04062
Phone: 1-800-377-9684
Fax: 1-207-892-2266
E-mail: Info@mainestandards.com