




TECHNICAL SERVICE BULLETIN
#10-003

Date Issued: June 01, 2010

Product:

VALIDATE [®] Product	REF	LOT	 <i>Extended Storage Temperature</i>
CM2	402bc	All Lots	-10° to -25°C
	402re	All Lots	-10° to -25°C
TDM1	301ab	All Lots	-10° to -25°C
	301bc	All Lots	-10° to -25°C
	301db	All Lots	-10° to -25°C
	301ri	All Lots	-10° to -25°C
	301vt	All Lots	-10° to -25°C
LP	501bc	All Lots	-10° to -25°C
	501db	All Lots	-10° to -25°C
	501ro	All Lots	-10° to -25°C
GC1	1100bc	All Lots	-10° to -25°C
	1100sd	All Lots	-10° to -25°C
	1100sa	All Lots	-10° to -25°C
GC2	1200bc	All Lots	-10° to -25°C
	1200sd	All Lots	-10° to -25°C
	1200sa	All Lots	-10° to -25°C
GC3	1300bc	All Lots	-10° to -25°C
	1300sd	All Lots	-10° to -25°C
	1300sa	All Lots	-10° to -25°C
GC4	1400bc	All Lots	-10° to -25°C
	1400sd	All Lots	-10° to -25°C
	1400sa	All Lots	-10° to -25°C

Subject: Extended Storage Temperature for VALIDATE[®] test kits

Explanation: Ongoing Quality Assurance testing has led us to extend the storage temperature for all lot numbers of VALIDATE[®] calibration verification / linearity test kits listed above from -10° to -20°C to -10° to -25°C.

Please retain this bulletin for your records.

Contact Information:

Customer Support Department
 Maine Standards Company
 765 Roosevelt Trail
 Windham, ME 04062
 Phone: 1-800-377-9684
 Fax: 1-207-892-2266
 E-mail: Info@mainestandards.com