





INSTRUCTIONS FOR COMPLETING THE DATA SUBMISSION FORM

In order for your data to be processed, all demographic information must be completed on at least one Data Submission Form submitted **per instrument**. A separate data submission form is required for each product and each instrument system.

Additional Data Submission Forms are available at www.mainestandards.com/analysis/dataformdown.php

If dilutions are performed, please use the Dilution Data Submission form also available at the address above or by calling Customer Service at 1-800-377-9684 or +1 207-892-1300.

Completed reports are typically returned within 5 business days from receipt of submitted data.

Step 1	1:	Record	l al	l necessary o	demograph	ic in	formatior	า in th	e spaces	provided.

Instrument system/model#: Acme Chemistry Analyzer Instrument ID/SN/Name SN 1049 Institution name: LGC Maine Standards Mailing address: 221 US Route 1 Cumberland Foreside, ME 04110 Attention Ruben Morse F-Mail: msc.datareduction@LGCGroup.com Contact name: **Data Reduction Department** Phone number: 1-800-377-9684 or 207-892-1300 Fax number: 207-892-2266

Step 2: Select how you would like to receive your completed calibration verification / linearity reports.

Select method of return:

[x] electronic copy via e-mail [] hardcopy via US Mail

Step 3: Record the testing date. Testing date: 03/21/18

Step 4: Record the product information with the required information.

Product: GC 1 *Set: Kit Lot #: 11AA08818 Kit Exp: 2019-10-15

*If using a kit with multiple sets

Step 5: Record the Method ID and reporting units for the analytes assayed. For information on determining the Method ID, please visit www.mainestandards.com/analysis/peer-group-setup-information.php.

Step 6: Record the duplicate or triplicate recovered values for analyte ID entered.

* Note: not all products contain Level 6.

A Method ID is REQUIRED to receive Peer Group Analysis.

	Analyte/Method ID:	BUN01	Units: mg/dL	
Level	Rep 1	Rep 2	Rep 3	
1	Ö	0	0	
2	38	37	38	
3	74	76	75	
4	112	113	113	
5	149	151	150	
6				

Step 7: Repeat Steps 4 thru 6 for all additional analytes assayed in the product.

Step 8: Repeat Steps 1 thru 7, on a separate Data Submission Form, for each additional VALIDATE® product to be submitted. Be sure to use a separate form for each instrument system.

Step 9: Submit completed Data Submission Forms to LGC Maine Standards for processing:

Email: msc.datareduction@LGCGroup.com (electronic forms)

Fax: 1-207-352-5253 US Mail: 221 US Route 1

Attn: Data Reduction Department Cumberland Foreside, ME 04110

If you need additional assistance completing the Data Submission Form, please call LGC Maine Standards Technical Support at 1-800-377-9684 +1 207-892-1300.