Product Announcement



March 2018

Dear Roche cobas[®] Customer,

Through an ongoing quality review process, LGC Maine Standards announces improvements to the VALIDATE[®] GC3re kit, order number 1300re. VALIDATE[®] GC3re has been reduced from a six level to a five level set. This change results in fewer samples to handle and fewer tests to run to validate the analytical measurement range (AMR) of the Roche cobas system enzymes.

The range coverage of VALIDATE[®] GC3re remains unchanged, despite the kit being reduced from six to five levels. Please see the chart below for the typical coverage of VALIDATE[®] GC3re.

Typical Values by Level 1300re						
Analyte	Units	1	2	3	4	5
ALP	U/L	5	304	603	901	1,200
ALT	U/L	4	178	352	526	700
AMY	U/L	3	502	1,002	1,501	2,000
AST	U/L	5	204	403	601	800
CK	U/L	7	505	1,004	1,502	2,000
GGT	U/L	3	302	602	901	1,200
LD	U/L	10	258	505	753	1,000
LIP	U/L	3	77	152	226	300

If you have any questions, please contact LGC Maine Standards by phone at either 1-800-377-9684 or +1-207-892-1300, or by e-mail at <u>msc.techsupport@LGCGroup.com</u>.

LGC Maine Standards deeply appreciates your business and we look forward to assisting your linearity / calibration verification needs in the future.

Sincerely,

Technical Support Department Maine Standards Company, LLC 221 US Route 1 Cumberland Foreside, ME 04110 Phone: 1-800-377-9684 ext. 1428 (United States and Canada) +1-207-892-1300 ext. 1428 Fax: +1-207-892-2266 E-mail: msc.info@LGCGroup.com