

PRODUCT ANNOUNCEMENT

<u>Date</u>: March 05, 2018 <u>Subject</u>: TDM1vt - Salicylate

Dear Customer,

Thank you for choosing VALIDATE® for your linearity and calibration verification needs. As a user of an Ortho Vitros analyzer, we would like to inform you of a product improvement to VALIDATE® TDM1vt (Order number: 301vt) which resulted from regular, ongoing quality reviews.

To improve salicylate's linear response on the Ortho Vitros systems, it has been removed from the original multi-analyte set and is now included as a single analyte in a second set of vials. The range coverage for salicylate and the remaining ten analytes in the first set all remain unchanged.

Please review the chart below for the typical coverage and the new configuration:

Typical Values for 301vt							
Set 1							
Analyte		Units	Level 1	Level 2	Level 3	Level 4	Level 5
Acetaminophen	ACTM	μg/mL	10	58	105	153	200
Carbamazepine	CARB	μg/mL	3.0	7.3	11.5	15.8	20.0
Digoxin	DIGN	ng/mL	0.40	1.3	2.2	3.1	4.0
Gentamicin	GENT	μg/mL	0.6	3.0	5.3	7.7	10.0
Phenobarbital	PHNO	μg/mL	3.0	22.3	41.5	60.8	80.0
Phenytoin	PHYT	μg/mL	3.0	12.3	21.5	30.8	40.0
Theophylline	THEO	μg/mL	1.0	10.8	20.5	30.3	40.0
Tobramycin	TOB	μg/mL	0.6	3.0	5.3	7.7	10.0
Valproic Acid	VALP	μg/mL	10	45	80	115	150
Vancomycin	VANC	μg/mL	5.0	16.3	27.5	38.8	50.0
Set 2							
Analyte		Units	Level 1	Level 2	Level 3	Level 4	Level 5
Salicylate	SALY	mg/dL	1.0	10.8	20.5	30.3	40.0

You do not need to make any changes – your next order of TDM1vt, order number: 301vt, will ship with this improved configuration. However, should you have any questions, please contact LGC Maine Standards Technical Support by phone at either 1-800-377-9684 (United States and Canada) or +1-207-892-1300 ext. 1428, or by e-mail at msc.techsupport@lgcgroup.com.

LGC Maine Standards continues to improve and expand our instrument-specific product offering for the Ortho Vitros, making linearity and calibration verification testing easier, with faster results laboratories can rely on. Click here to view our complete offering for your analyzer (which now includes Vitamin D - 506vt).

We appreciate your business and look forward to assisting you in the future.

Sincerely,

Technical Support Department LGC Maine Standards msc.techsupport@lgcgroup.com

MAINE STANDARDS COMPANY, LLC
221 US ROUTE 1 • CUMBERLAND FORESIDE, MAINE 04110

TEL: 207.892.1300 • FAX: 207.892.2266 • TOLL-FREE: 800.377.9684

www.mainestandards.com