

TECHNICAL SERVICE BULLETIN

2019-011



Date Issued: March 22, 2019

Subject: Removal of 6th Level in GC1 (1100ro)

Explanation:

LGC Maine Standards has removed the 6th level in GC1 (1100ro). The 6th level was included in the GC1 Kit to meet the analytical measuring range for the following analytes: Albumin, BUN and Triglycerides, on the Roche Hitachi Mod P and D analyzers only. Since Roche has discontinued support of the Hitachi Mod P and Mod D instruments (see Roche Diagnostics Announcement 17-184, dated 8-8-17), the 6th level is no longer needed. GC1 (1100ro) Levels 1 through 5 will continue to meet the analytical measuring range of all included analytes for the Roche cobas[®] systems.

Please contact LGC Maine Standards Technical Support if you have any questions at 1-800-377-9684.

Contact Information: Technical Support Department
Maine Standards Company, LLC
221 US Route 1
Cumberland Foreside, ME 04110
Phone: 1-800-377-9684 ext. 1428
+1-207-892-1300 ext. 1428
Fax: +1-207-892-2266
E-mail: msc.technicalsupport@LGCGroup.com