TECHNICAL SERVICE BULLETIN 19-001





Date Issued: March 14, 2019

Subject: LGC Maine Standards' Shipping Protocol

Explanation: LGC Maine Standards products shipped with frozen gel packs and indicated

with "Cool Arrival" on the packing slip, may arrive thawed when received at your facility. Thawing of product during shipment will not adversely impact product performance. Once product arrives at your facility, immediately store

at the recommended storage temperature.

VALIDATE® Heparin and VALIDATE® Fibrinogen products should arrive at your facility frozen and immediately be stored at the recommended storage

temperature.

Contact Information: Technical Support Department

Maine Standards Company, LLC

221 US Route 1

Cumberland Foreside, ME 04110 Phone: 1-800-377-9684 ext. 1428

+1-207-892-1300 ext. 1428

Fax: +1-207-892-2266

E-mail: msc.techsupport@LGCGroup.com