

TECHNICAL SERVICE BULLETIN 2021-002

Date Issued: February 3, 2021

Subject:

VALIDATE[®] GC2 1200vt Kit Labels Incorrectly List Uric Acid (UA)

VALIDATE®	REF	LOT
VALIDATE [®] GC2	1200vt	12AK095200

Explanation:

Customers who have received VALIDATE[®] GC2 1200vt lot 12AK095200 prior to January 28, 2021 have an incorrect label on the kit box which states that CO_2 , ETOH, NH₃, and UA are included in the product. Uric Acid (UA) is not included in the product. The analytes within the 1200vt product are CO_2 , ETOH, and NH₃.

We apologize for any inconvenience this may have caused your laboratory.

Contact Information:	Technical Support Department Maine Standards Company, LLC 221 US Route 1 Cumberland Foreside, ME 04110	
	Fax:	1-800-377-9684 ext. 1428 +1-207-892-1300 ext. 1428 +1-207-892-2266 msc.techsupport@LGCGroup.com