TECHNICAL SERVICE BULLETIN 21-006



Date Issued: May 20, 2021

Subject: Gentamicin in 301ri Kit Configuration

VALIDATE® TDM1 (Order Number 301ri) test kit is currently being manufactured using the GentQ configuration on the Roche cobas®.

The Roche cobas® currently has two Gentamicin configurations available, the GentQ and the newer version, GentC. At this time the current VALIDATE® TDM1 (301ri) kit will recover half of the new GentC reagent reportable range (0.24-12.0 ug/mL) due to the reportable range and method difference between GentQ and GentC. The GentQ reportable range is 0.3-10.0 ug/mL.

To address the GentC range coverage, we have a set 2 in development. The set 2 for the VALIDATE® TDM1, 301ri, kit will be specific for the GentC reagent and associated reportable range. We anticipate this kit to be available in September 2021.

For questions or additional information, please contact the LGG Maine Standards Technical Support department.

Contact Information: Technical Support Department

Maine Standards Company, LLC

221 US Route 1

Cumberland Foreside, ME 04110 Phone: 1-800-377-9684 ext. 1428 +1-207-892-1300 ext. 1428

11-207-092-1000 6

Fax: +1-207-892-2266

E-mail: <u>msc.info@LGCGroup.com</u>