

TECHNICAL SERVICE BULLETIN 2024-018



Date Issued: January 27, 2025

Subject: Safety Data Updates for VALIDATE Products

Explanation: Safety Data Sheets (SDS), Labeling and Product Inserts for VALIDATE Products have been updated to reflect the latest EU, OSHA, and GHS requirements.

VALIDATE TM (407RE) and (407AB) SDS, Labels and Product Insert has been updated for information related to Triton X-100. Triton X-100 is on the European Union and United Kingdom's list for substances of very high concern (SVHC) which must be identified in SDS and package inserts if present above 0.1% in products/materials supplied to the UK or EU. Please see SDS-000044 for additional information.

VALIDATE UC4 (704xx) and UC6 (706af) SDS, Labels and Product Inserts have been updated for information related to corrosive materials. The following symbols have been added to SDS-000002 and SDS-000034:



VALIDATE ACTH (302MM), SDS-000045 has been updated to include a GHS harmful to health and exclamation hazard:



The following Warning symbol has been added to Outer Box Labeling for the products above:



Contact Information:

Technical Support Department
LGC Clinical Diagnostics, Inc.
221 US Route 1
Cumberland Foreside, ME 04110
Phone: 1-800-377-9684
+1-207-892-1300
Fax: +1-207-892-2266
E-mail: msc.techsupport@LGCGroup.com