

# TECHNICAL SERVICE BULLETIN

## 2024-016



**Date Issued:** Jun 24, 2024

**Subject:** VALIDATE® GC 4 (1400re) Lot: 10706723 forward

**Explanation:** Updated recovery targets for DBIL

VALIDATE® GC4 (1400re), for Roche cobas®, **now targets the DBIL2 method's measuring range.** Updated targets for DBIL2 will be implemented on lot number 10706723 and will be used going forward. Please update your records and/or staff appropriately as recoveries will vary compared to previous lots.

Please see the chart below for the typical coverage of VALIDATE® GC4 (1400re),

Analyte	Units	Analyzer Range	Validate Range	Percent Coverage
DBIL	mg/mL	0.08 – 13.8	0.1 – 15.5	0.1% – 112%
TBIL	mg/mL	0.15 – 35.1	0.2 – 37.5	0.1% - 107%

**Contact Information:** Technical Support Department  
LGC Clinical Diagnostics, Inc.  
221 US Route 1  
Cumberland Foreside, ME 04110  
Phone: 1-800-377-9684  
+1-207-892-1300  
Fax: +1-207-892-2266  
E-mail: [msc.techsupport@LGCGroup.com](mailto:msc.techsupport@LGCGroup.com)