

TECHNICAL SERVICE BULLETIN 2025-005



Date Issued: Mar 25, 2025

Subject: VALIDATE® UC4 (704ab) Method Updates

Explanation:

We are pleased to announce that VALIDATE® UC4 (704ab) is now targeting Abbott's new Crea2 and Phos2 assays.

Please see the chart below for the typical coverage of VALIDATE® UC4 (704ab).

Analyte	Units	Analyzer Range	VALIDATE Range	Percent Coverage
CA	mg/dL	2 – 24	3.3 – 22	6% – 91%
CREA	mg/dL	2.56 – 740	7 – 667	0.6% – 90%
MG	mg/dL	1.81 – 26.35	3.5 – 24	7% – 90%
PHOS	mg/dL	2.4 – 186.2	5 – 117	1% – 95%

Please contact LGC Maine Standards Technical Support if you have any questions at 1-800-377-9684.

Contact Information: Technical Support Department
LGC Clinical Diagnostics, Inc.
221 US Route 1
Cumberland Foreside, ME 04110
Phone: 1-800-377-9684
+1-207-892-1300
Fax: +1-207-892-2266
E-mail: msc.techsupport@LGCGroup.com