

TECHNICAL SERVICE BULLETIN 25-008



Date Issued: Apr 11, 2025

Subject: VALIDATE® GC4 (1400ab) Method Update

Explanation: We are pleased to announce that VALIDATE® GC4 (1400ab) is now targeting Abbott's new TBili2 assay.

Please see the chart below for the typical coverage of VALIDATE® GC4 (1400ab).

| Analyte | Units | Analyzer Range | VALIDATE Range | Percent Coverage |
|---------|-------|----------------|----------------|------------------|
| DBIL | mg/dL | 0.1 – 15 | 0.1 – 14.1 | 0% – 94% |
| TBIL | mg/dL | 0.1 – 25 | 0.2 – 23.7 | 0.4% – 95% |

Please contact LGC Maine Standards Technical Support if you have any questions at 1-800-377-9684.

Contact Information: Technical Support Department
LGC Clinical Diagnostics, Inc.
221 US Route 1
Cumberland Foreside, ME 04110
Phone: 1-800-377-9684
+1-207-892-1300
Fax: +1-207-892-2266
E-mail: msc.techsupport@LGCGroup.com