

# TECHNICAL SERVICE BULLETIN 25-011



**Date Issued:** Apr 16, 2025

**Subject:** VALIDATE<sup>®</sup> CM2 (402vt) Design change

**Explanation:** We are pleased to announce that VALIDATE<sup>®</sup> CM2 (402vt) is now targeting the Sentinel Diagnostics CRP Ultra (MP) assay using the Standard method.

- Standard Method users: Use Levels 1-5.
- Ultrasensitive Method users: Use Levels 1-3.

Please see the chart below for the typical coverage of VALIDATE<sup>®</sup> CM2 (402vt).

Analyte	Units	Analyzer Range	VALIDATE Range	Percent Coverage
hsCRP	mg/L	0.3 – 320	0.5 – 284	0.1% – 89%
NTproBNP	pg/mL	20 – 30000	29 – 25703	0.03% – 86%
Tnl	ng/mL	0.012 – 80	0.05 – 70	0.05% – 87%

Please contact LGC Maine Standards Technical Support if you have any questions at 1-800-377-9684.

**Contact Information:** Technical Support Department  
LGC Clinical Diagnostics, Inc.  
221 US Route 1  
Cumberland Foreside, ME 04110  
Phone: 1-800-377-9684  
+1-207-892-1300  
Fax: +1-207-892-2266  
E-mail: [msc.techsupport@LGCGroup.com](mailto:msc.techsupport@LGCGroup.com)