

TECHNICAL SERVICE BULLETIN

25-014



Date Issued: Jun 9, 2025

Subject: AccuTrak™ Calibration Verification Panel Announcement

Explanation:

We are pleased to announce that LGC Maine Standards has developed a new AccuTrak™ Calibration Verification Panel product line.

AccuTrak™ vs. VALIDATE®

AccuTrak™ Calibration Verification Panels are designed as individual levels to provide coverage across the reportable range of your assay (please see product list instrument specific designs and best-fits) and meet Calibration Verification requirements found in CLIA and ISO 15189. Products are liquid, ready-to-use, good for multiple uses with extended open-vial stability and leverage MSDRx® Infinity for peer data comparison.

VALIDATE® Calibration Verification and Linearity products (Levels 1-5) are designed according to CLSI's EP06 guidelines in an equal-delta formulation and meet Calibration Verification and Linearity requirements outlined in CLIA and ISO 15189. Products are liquid, ready-to-use, good for multiple uses with extended open-vial stability and leverage MSDRx® Infinity for peer data comparison and linearity analysis

AccuTrak™	VALIDATE®
For calibration verification and reportable range verification	For linearity, calibration verification, and reportable range verification
Single analyte	Single or Multi-analyte
Liquid, ready-to-use	Liquid, ready-to-use
Good for multiple uses and have extended open-vial stability	Good for multiple uses and have extended open-vial stability
Configurations tailor-made to analytical instrument platforms	Configurations tailor-made to analytical instrument platforms
Levels are prepared individually across the reportable range	Levels 1-5 are prepared according to CLSI's EP06-A guidelines
Data reduction software and peer group comparison through MSDRx Infinity	Data reduction software and peer group comparison through MSDRx Infinity

Please contact LGC Maine Standards Technical Support if you have any questions at 1-800-377-9684.

Contact Information:

Technical Support Department
LGC Clinical Diagnostics, Inc.
221 US Route 1
Cumberland Foreside, ME 04110
Phone: 1-800-377-9684
+1-207-892-1300
Fax: +1-207-892-2266
E-mail: msc.techsupport@LGCGroup.com